

Omniwatch OW410 Nurse Call Monitor



Provides call light reporting at lowest cost. Also can send an email to a cell phone or speak over a two way radio with the call information.

Wires across the annunciator lamp of existing call systems. Requires 12 - 24 volts AC or DC from the call system.

This unit will read the data from the call system and send it to a computer. The computer will display the active calls. The computer has a basic call reporting function. There is also an option to use the full Omniwatch reporting program. The computer can also send an email to a cell phone or speak over a two way radio with the room that placed the call.

A computer needs to be dedicated to the nurse call system. You can use one of your computers or buy one from Omniwatch.

The OW410 connects to the computer with a RS232 or RS422 serial cable. There is also an option to use WIFI wireless networking for the connection.

The software reporting program will display a list of calls with response time for a given date range. The list of calls can also be for only a single room. The total of calls in the list and the average response time is also displayed.